



# KALOYA MICRONEEDLING PEN

Advanced Skin Rejuvenation System

*Kaloya*<sup>™</sup>  
DERMATECH



Health  
Canada

Santé  
Canada


# Key Features



Adjustable needle  
depth: 0.25 mm –  
2.0 mm



6 speed levels up  
to 7700 RPM



Cordless,  
lightweight (93g)  
design



Rechargeable with  
LED charge  
indicator



Single-use EO  
sterilized  
cartridges



Protective sleeve  
for hygiene

# Solution for These Skin Conditions

Acne scars


Fine lines and wrinkles

Enlarged pores

Uneven skin tone and texture

Stretch marks and mild pigmentation


Surgical scars



# Before & After Results




# How to Use



# Depth Guidelines

ORBITAL RIM: 0.25–0.5 MM (AVOID DIRECT CONTACT)



NOSE: 0.25–0.75 MM



FOREHEAD: 0.25–1.0 MM



ACNE SCARS: UP TO 1.5 - 2.5MM



# Certified Training Included

Hands-on  
microneedling  
training with every  
purchase

Covers theory,  
safety protocols,  
and depth  
selection

Includes treatment  
planning and  
contraindications

Certified by Kaloya  
Dermatech Ltd

Delivered by  
experienced  
professionals

---

# Safety & Compliance Standards


---

- Complies with ISO 10993-1, IEC60601-1, and IEC62366 standards
- EO sterilized, single-use cartridges
- Safe for all Fitzpatrick skin types I–VI
- Adjustable depth for precision treatments
- Protective sleeve prevents cross-contamination



# Ideal For Professionals Who:

- Medspa, skincare clinic, and dermatology professionals
- Advanced facial and anti-aging service providers
- Businesses looking to expand treatment menus
- Professionals seeking high-quality, certified equipment
- Mobile service providers needing portable tools





---

# Contact & Ordering

-  Email: [info@kaloyadermatech.com](mailto:info@kaloyadermatech.com)
-  Phone: (604) 771-6450
-  Website: [www.kaloyadermatech.com](http://www.kaloyadermatech.com)

**Book a demo or place an order today!**